

# Human IgMλ Isotype Control

Cat. No.	Format	Size
0158L-01	Purified (UNLB)	0.5 mg

#### **Description**

**Isotype** Human IgMλ

Source Human myeloma serum

### **Applications**

Quality tested applications include -

**ELISA** 

## **Working Dilutions**

**ELISA** Purified (UNLB) antibody  $\leq 2 \mu g/mL$ 

Other Applications Since applications vary, you should determine the optimum working dilution for the product that is

appropriate for your specific need.

#### Handling and Storage

- The purified (UNLB) IgM is provided at a concentration of 0.5 mg/mL in 1.0 mL of borate buffered saline, pH 8.2. No preservatives or amine-containing buffer salts added. Store at 2-8°C.
- Reagents are stable for the period shown on the label if stored as directed.
- If turbidity occurs, clarify the solution by centrifugation prior to use.

In accordance with current Good Manufacturing and Good Laboratory Practices (cGMP/cGLP), any protein of human blood origin should be handled pursuant to your organization's documented safety procedures and as if it is capable of transmitting infection. This product has **NOT** been tested for viral, bacterial, or other infectious agents such as, but not limited to, HIV, HbsAq, and HCV.

For Research Use Only. Not for Diagnostic or Therapeutic Use.

TB0158L 24-Mar-15

Email: info@southernbiotech.com • Website: www.southernbiotech.com