



## TMB Stop Solution (650 nm)

Cat. No.	Format	Size
0413-01	Solution	100 mL
0413-01L	Solution	1 L

### Overview

---

**Reaction Product** Stops the reaction of TMB One Component Microwell Substrate, Solution (SB Cat. No. 0410-01 & 0410-01L) without a spectral change from blue (650 nm) to yellow (450 nm)

### Applications

---

ELISA

### Product Use

---

- Ready to use formulation
- Use equal volume of TMB substrate and TMB Stop Solution

### Handling and Storage

---

- Recommended storage is at 25°C. Refrigerated temperatures will not harm the reagent.
- Reagent is stable for the period shown on the label if stored as directed.

***For Research Use Only. Not for Diagnostic or Therapeutic Use.***

TB0413  
24-Jul-19

---

**Corporate Offices:** 160 Oxmoor Blvd • Birmingham, AL 35209 • USA **Mailing Address:** P.O. Box 26221 • Birmingham, AL 35260 • USA

**Tel:** 205.945.1774 • U.S. and Canada: 800.722.2255 • **Fax:** 205.945.8768

**Email:** [info@southernbiotech.com](mailto:info@southernbiotech.com) • **Website:** [www.southernbiotech.com](http://www.southernbiotech.com)