
SouthernBiotech =



Goat Anti-Mouse IgA

Cat. No.	Format	Size
1040-01	Purified (UNLB)	1.0 mg
1040-02	Fluorescein (FITC)	1.0 mg
1040-03	Rhodamine (TRITC)	1.0 mg
1040-04	Alkaline Phosphatase (AP)	1.0 mL
1040-05	Horseradish Peroxidase (HRP)	1.0 mL
1040-07	Texas Red® (TXRD)	1.0 mg
1040-08	Biotin (BIOT)	1.0 mg
1040-09	R-phycoerythrin (PE)	0.5 mg
1040-30	Alexa Fluor® 488 (AF488)	1.0 mg
1040-31	Alexa Fluor® 647 (AF647)	1.0 mg
1040-32	Alexa Fluor® 555 (AF555)	1.0 mg



FLISA plate was coated with purified mouse IgA, IgG, and IgM. Immunoglobulins were detected with serially diluted Goat Anti-Mouse IgA-PE (SB Cat. No. 1040-09).

Description

Specificity Reacts with the heavy chain of mouse IgA

Pooled antisera from goats hyperimmunized with mouse IgA Source

Cross Adsorption Mouse IgG and IgM; may react with immunoglobulins from other species **Purification** Affinity chromatography on mouse IgA covalently linked to agarose

Applications

Quality tested applications include -

ELISA 1-15,23

FLISA FC ^{1,5,6,10,13,17-20,23,24}

Other referenced applications include -

ELISPOT 2,7,14,16,23 IHC-FS 3,4,9,11 IHC-PS 5,6 ICC 21,24 WB 8,12,22

Sep ¹

Depletion ²³

Working Dilutions

ELISA	AP conjugate	1:2,000 – 1:4,000
	HRP conjugate	1:4,000 - 1:8,000
	BIOT conjugate	1:5,000 - 1:20,000
FLISA	FITC, TRITC, TXRD, AF488, and AF555 conjugates	1:100 – 1:400
	PE and AF647 conjugates	≤ 1 μg/mL
Flow Cytometry	FITC, BIOT, and AF488 conjugates	\leq 1 μ g/10 ⁶ cells
	PE and AF647 conjugates	≤ 0.1 μg/10 ⁶ cells
	For flow cytometry, the suggested use of these reagents is in a final volume of 100 μ L	
Other Applications	Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.	

For Research Use Only. Not for Diagnostic or Therapeutic Use.

Handling and Storage

- The purified (UNLB) antibody is supplied as 1.0 mg purified immunoglobulin in 1.0 mL of borate buffered saline, pH 8.2. No preservatives or amine-containing buffer salts added. Store at 2-8°C.
- The fluorescein (FITC), rhodamine (TRITC), and Texas Red® (TXRD) conjugates are supplied as 1.0 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The alkaline phosphatase (AP) conjugate is supplied as 1.0 mL in a stock solution of 50 mM Tris/1 mM MgCl₂/50% glycerol, pH 8.0, containing NaN₃ as preservative. Store at 2-8°C or long-term at -20°C.
- The horseradish peroxidase (HRP) conjugate is supplied as 1.0 mL in a stock solution of 50% glycerol/50% PBS, pH 7.4. No preservative added. Store at 2-8°C or long-term at -20°C.
- The biotin (BIOT) conjugate is supplied as 1.0 mg in 2.0 mL of PBS/NaN₃. Store at 2-8°C.
- The Alexa Fluor® 488 (AF488), Alexa Fluor® 555 (AF555), and Alexa Fluor® 647 (AF647) conjugates are supplied as 1.0 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The R-phycoerythrin (PE) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN₃ and a stabilizing agent. Store at 2-8°C. **Do not** freezel
- Protect fluorochrome-conjugated forms from light. Reagents are stable for the period shown on the label if stored as directed.

Warning

Some reagents contain sodium azide. Please refer to product specific SDS.

References

- Harriman GR, Kunimoto DY, Elliott JF, Paetkau V, Strober W. The role of IL-5 in IgA B cell differentiation. J Immunol. 1988;140:3033-9. (ELISA, FC, Sep)
- van Ginkel FW, Wahl SM, Kearney JF, Kweon M, Fujihashi K, Burrows PD, et al. Partial IgA-deficiency with increased Th2-type cytokines in TGF-β1 knockout mice. J Immunol. 1999;163:1951-7. (ELISA, ELISPOT) 2.
- Pappo J, Torrey D, Castriotta L, Savinainen A, Kabok Z, Ibraghimov A. Helicobacter pylori infection in immunized mice lacking major histocompatibility complex class I and class II functions. Infect Immun. 1999;67:337-41. (ELISA, IHC-FS) 3.
- Dunn LA, Upcroft JA, Fowler EV, Matthews BS, Upcroft P. Orally administered Giardia duodenalis extracts enhance an antigen-specific antibody response. Infect Immun. 2001;69:6503-10. (ELISA, IHC-FS)
- Ibnou-Zekri N, Blum S, Schiffrin EJ, von der Weid T. Divergent patterns of colonization and immune response elicited from two intestinal Lactobacillus strains that display similar properties in vitro. Infect Immun. 2003;71:428-36. (ELISA, FC, IHC-PS)
- Wang J, Anders RA, Wu Q, Peng D, Cho JH, Sun Y, et al. Dysregulated LIGHT expression on T cells mediates intestinal inflammation and contributes to IgA nephropathy. J Clin Invest. 2004;113:826-35. (ELISA, FC, IHC-PS)
- Hovden A, Brokstad KA, Major D, Wood J, Haaheim LR, Cox RJ. A pilot study of the immune response to whole inactivated avian influenza H7N1 virus vaccine in mice. Influenza Other Respir Viruses. 2009;3:21-8. (ELISA, ELISPOT)
- Sanchez CJ, Kumar N, Lizcano A, Shivshankar P, Hotopp JC, Jorgensen JH, et al. Streptococcus pneumoniae in biofilms are unable to cause invasive disease due to altered virulence 8. determinant production. PLoS One. 2011;6(12):e28738. (ELISA, WB)
- Harris DP, Vogel P, Wims M, Moberg K, Humphries J, Jhaver KG, et al. Requirement for class II phosphoinositide 3-kinase C2α in maintenance of glomerular structure and function. Mol Cell Biol. 2011;31:63-80. (ELISA, IHC-FS)
- Willems KN. Regulation of humoral immunity by Pim kinases [dissertation]. Worcester (USA): University of Massachusetts; 2011. (ELISA, FC)
- Zhang H, Jia C, Xi H, Li S, Yang L, Wang Y. Specific inhibition of Candida albicans growth in vitro by antibodies from experimental Candida keratitis mice. Exp Eye Res. 2011;93:50-8. (ELISA, IHC-FS)
- 12 Nakanishi K, Narimatsu S, Ichikawa S, Tobisawa Y, Kurohane K, Niwa Y, et al. Production of hybrid-lqG/lgA plantibodies with neutralizing activity against Shiga toxin 1. PLoS One. 2013;8(11):e80712. (ELISA, WB)
- Goverse G, Olivier BJ, Molenaar R, Knippenberg M, Greuter M, Konijn T, et al. Vitamin A metabolism and mucosal immune function are distinct between BALB/c and C57BL/6 mice. Eur J Immunol. 2015;45:89-100. (ELISA, FC)
 Surman SL, Jones BG, Rudraraju R, Sealy RE, Hurwitz JL. Intranasal administration of retinyl palmitate with a respiratory virus vaccine corrects impaired mucosal IgA response in
- the vitamin A-deficient host. Clin Vaccine Immunol. 2014; 21:598-601. (ELISA, ELISPOT) Pucella JN, Yen W, Kim MV, van der Veeken J, Socci ND, Naito Y, et al. miR-182 is largely dispensable for adaptive immunity: lack of correlation between expression and function. J
- Immunol. 2015;194:2635-42. (ELISA) Patel PS, Kearney JF. Neonatal exposure to pneumococcal phosphorylcholine modulates the development of house dust mite allergy during adult life. J Immunol. 2015;194:5838-50. (FLISPOT)
- 17. Jiang H, Bos NA, Cebra JJ. Timing, localization, and persistence of colonization by segmented filamentous bacteria in the neonatal mouse gut depend on immune status of mothers and pups. Infect Immun. 2001;69:3611-7. (FC)
- Cook AJ, Oganesian L, Harumal P, Basten A, Brink R, Jolly CJ. Reduced switching in SCID B cells is associated with altered somatic mutation of recombined S regions. J Immunol. 2003;171:6556-64. (FC)
- Guikema JE, Schrader CE, Leus NG, Ucher A, Linehan EK, Werling U, et al. Reassessment of the role of Mut S homolog 5 in Ig class switch recombination shows lack of involvement in cis- and trans-switching. J Immunol. 2008;181:8450-9. (FC)
- Pauklin S, Petersen-Mahrt SK. Progesterone inhibits activation-induced deaminase by binding to the promoter. J Immunol. 2009;183:1238-44. (FC)
- Armitage CW, O'Meara CP, Harvie MC, Timms P, Wijburg OL, Beagley KW. Evaluation of intra- and extra-epithelial secretory IgA in chlamydial infections. Immunology. 2014;143:520-
- Endt K, Stecher B, Chaffron S, Slack E, Tchitchek N, Benecke A, et al. The microbiota mediates pathogen clearance from the gut lumen after non-typhoidal Salmonella diarrhea. PLoS Pathog. 2010;6(9):e1001097. (WB)
 Harriman GR, Kunimoto DY, Elliott JF, Paetkau V, Strober W. The role of IL-5 in IgA B cell differentiation. J Immunol. 1988;140:3033-9. (ELISA, FC, ELISPOT, Depletion)
- 23.
- Usui T, Wakatsuki Y, Matsunaga Y, Kaneko S, Koseki H, Kita T, et al. Overexpression of B cell-specific activator protein (BSAP/Pax-5) in a late B cell is sufficient to suppress differentiation to an Ig high producer cell with plasma cell phenotype. J Immunol. 1997;158:3197-204. (FC, ICC)

Texas Red® is a registered trademark of Molecular Probes, Inc.

Alexa Fluor[®] 488, 847, and 555 are provided under an Intellectual property license from Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product purchased annotation to the product and components or the product or its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. USA or outlicensing@lifetech.com. For information on purchasing a license to this product for any other use, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008

> TB1040 05-Mar-18